

# AUDITEST

## SUPER PORTABLE WITH A PROFESSIONAL CLASS



The AudiTest, screening manual audiometer, has been designed for achieving an accurate diagnosis of hearing faculties, wherever you are. Modern and performant, its size and ergonomics allow the health's professional to use it easily with just one hand.



- 11 frequencies, from 125 Hz up to 8000 Hz
  - Sound levels adjustable from -10 up to 100 HLdB with 5 dB steps
  - Continuous or pulsed tone
  - Fedded with 2 single 1,5 V battery sticks
  - Display on a large LCD screen :
    - Frequency and level
    - Stimulus mode
    - Status of batteries
    - Sound delivery indicator
  - Silent and embossed keys
  - Automatic power shut off
  - Medical CE marking (G-MED / CE 0459)
- In compliance with audiometry standard EN 60645-1

- Delivered into a kit, including :
- Audiometric box
  - Headphones
  - 2 batteries
  - Instructions for use and maintenance

**WARRANTY : 3 YEARS**  
**FRENCH DESIGN AND MANUFACTURING**

Type	4 – screening (meets and exceeds EN 60645-1)
Output	CA – Air Conduction
Frequency range	125, 250, 500, 750, 1000, 1500, 2000, 3000, 4000, 6000, 8000 Hz
Output level range	- 10 up to 100 HLdB
Sound level pitch	5 dB
Stimulus	Continuous, pulsed
Headphones	Sennheiser: light, comfortable, highly isolated
Power supply	2 batteries of 1,5 V – type AA / LR06
Accuracy	Frequencies +/- 2 % - Levels +/- 3 dB
Audiometric box dimensions	79 x 122 x 33 mm
Total weight of the kit	600 g (155 g for the audiometric box)